FebriDx[®]

| | | |
|----|------------------|----|
| EN | English | 2 |
| ES | Español | 7 |
| NL | Nederlands | 13 |
| FR | Français | 19 |
| DE | Deutsche | 25 |
| IT | Italiano | 31 |

INTENDED USE

FebriDx is a rapid immunoassay for the visual, qualitative, in vitro detection of elevated levels of both MxA and CRP directly from fingerstick whole blood. The test measures a clinically significant immune response to a suspected invasive viral and/or bacterial infection in patients older than 2 years that present within 3 days of an acute onset fever (exhibited or reported) and within 7 days of new onset respiratory symptoms consistent with a community-acquired upper respiratory infection.

The FebriDx test aids in the clinical identification of patients with an underlying invasive viral infection from either Influenza A/B, Adenovirus, Respiratory Syncytial Virus, Metapneumovirus, Parainfluenza Virus, or Epstein-Barr Virus; and/or patients with a clinically significant immune response consistent with an underlying bacterial infection.

The test is intended for professional use in an outpatient setting and should be used in conjunction with other clinical evidence including laboratory, radiographic, and epidemiological information.

Negative results do not preclude respiratory infection (e.g. rhinovirus, coronavirus) and should not be used as the sole basis for diagnosis, treatment, or other clinical and patient management decisions. In addition to utilizing radiography and clinical presentation to aid in diagnosis, additional laboratory testing (e.g. bacterial and viral culture, immunofluorescence, and polymerase chain reaction [PCR]) may be used to confirm whether a specific respiratory pathogen exists.

Acute respiratory infections (ARIs) affect 20% of the population annually. The significant overlap in symptoms and signs makes it challenging for physicians to differentiate viral from bacterial infection and to identify which patients require antibiotic therapy. More than 50% of antibiotics prescribed for ARIs are unnecessary as the vast majority are caused by viruses, in which antibiotics provide no clinical benefit to the patient.

BIOMARKERS

MxA

MxA becomes elevated in the presence of acute viral infection, but is not specific to a particular type of virus. MxA has a low basal concentration of less than 15 ng/mL, a fast induction time of 1-2 hours, and a long half-life of 2.3 days, making it an ideal marker for viral infection.¹ Numerous clinical studies demonstrate that MxA protein expression in peripheral blood has been shown to be a sensitive and specific marker for viral infection.^{1,6}

CRP

C-reactive protein (CRP) is a nonspecific, acute-phase protein that increases during an inflammatory process, especially following severe infection. Bacterial infection is a potent stimulus of marked CRP elevation, which occurs within 4-6 hours of infection and peaks after 36 hours.^{7,8} Some viral infections may cause CRP to elevate.^{5,10} The administration of antibiotics causes CRP levels to fall rapidly.^{11,12}

Multiplexed Pattern of Results

In isolation, neither MxA nor CRP alone is sensitive or specific enough to differentiate viral from bacterial infection. At low levels, CRP is very sensitive but non-specific at confirming a bacterial infection. At high levels, CRP becomes very specific for bacterial infection but has low sensitivity. MxA is specific for viral infection only and is insensitive for the presence of a bacterial infection. The FebriDx test produces a multiplexed pattern of results by simultaneously detecting elevated levels of MxA and CRP together to help identify patients suffering from clinically significant acute respiratory infection as well as differentiate viral from bacterial infectious etiology.

PRINCIPLES OF THE TEST

The FebriDx test is a lateral flow immunoassay. The test strip, which is located inside an all-in-one plastic housing that incorporates a built-in retractable lancet, blood collection and transfer tube, and button release buffer activation mechanism, utilizes monoclonal anti-MxA and anti-CRP antibodies to create a control line and two result lines (MxA and CRP). The cut-off value for MxA is 40 ng/mL and the cut-off value for CRP is 20 mg/L serum equivalent.

If the fingerstick blood samples contain elevated levels of MxA or CRP above their respective cut-off levels, the appropriate test line will appear in the Result Window. FebriDx is a disposable, rapid test requiring 10 minutes for a result.

Materials Provided

- 25 single use tests
- 1 package insert

Materials Not Provided

- Timer
- Alcohol
- Gauze
- Sterile dressing

WARNINGS AND PRECAUTIONS



1. For in vitro diagnostic use only.
2. Keep the FebriDx test in the sealed foil pouch until just before use.
3. Do not use the FebriDx test past the expiration date.
4. Use standard precautions for collecting and handling a blood sample.
5. All specimens should be considered potentially hazardous and handled in the same manner as an infectious agent.

6. Wear disposable gloves while handling samples and wash hands after the test is complete.
7. The lancet is sterile until the protective tab is removed. Do not use the lancet if the protective tab is not secured in place.
8. The FebriDx test is designed to proceed in sequential order and locking mechanisms exist to prevent skipping the prior step.
9. The FebriDx test is a single-use item with no reusable components. Proper handling and disposal methods should be established according to local, state, and federal regulations.
10. The FebriDx test requires a visual readout. Do not interpret the test result if you have color-impaired vision.
11. A brightly lit environment is recommended for interpreting the test results.


STORAGE AND STABILITY

Store the FebriDx test between 77°F/25°C and 39°F/4°C. The FebriDx tests are stable until the expiration dates printed on their packaging.


TEST COMPONENTS


Lance Finger

Collect Blood

Deliver Blood

Deliver Buffer

TEST PROCEDURE - Check the expiration date on all packaging.

1. Tear open the foil pouch at the indicated perforation and remove the test.

Collecting and Transferring the Fingerstick Blood Sample

Note: Use standard precautions for collecting and handling a blood sample.


2. Cleanse the fingertip with an alcohol pad and allow it to air dry.
3. Locate the lancet and remove the **Protective Lancet Tab**. Firmly press the lancet to puncture the skin. Wipe away the first drop of blood with gauze and gently massage towards the puncture site to encourage blood flow.
4. Place the **Blood Collection Tube** at approximately a 45 degrees angle below the finger blood sample as contact is made. Fill the **Blood Collection Tube** in its entirety by touching the blood to the tip of the **Blood Collection Tube**. If the **Blood Collection Tube** is not full, gently squeeze the finger and add more blood.

Note: Capillary action will automatically draw the blood sample into the Blood Collection Tube in the required amount (5 µl).

5. Once the **Blood Collection Tube** is filled with blood, rotate it over the **Blood Transfer Zone** to deliver it to the test. The **Blood Collection Tube** will lock into position. Wait for most of the blood to be transferred to test strip (~5-10 seconds) before proceeding to step #6 and activating the test.

Note: If the blood does not immediately begin to transfer the blood onto the test strip, reverse the **Blood Collection Tube's** rotation back to its original position. This will occur against some resistance. Add additional fingerstick blood to the **Blood Collection Tube** to ensure the tube is completely filled before rotating back onto the test strip.

6. Activate the test by firmly and fully pressing the **Buffer Release Button** to deliver the buffer. The **Buffer Release Button** should be pressed within 2 minutes of transferring the blood sample. If no fluid is visible within 25–30 seconds, firmly re-press the **Buffer Release Button**. Lay the test on a flat surface and wait for 10 minutes. Results will appear in the **Result Window**.

TEST RESULTS

Results should be interpreted at the 10 minute mark. Do not read the test results after 1 hour.

An unused test, or test that has not yet been activated by pressing the **Buffer Release Button**, will show three faint orange lines in the Result Window.

A blue control line must appear in the **Result Window** for the test to be valid.

Positive Result

The positive result lines appear as red or black lines in the **Result Window**. An uneven or incomplete result line is due to an uneven sample distribution on the test strip. Even if the result line is faint in color, incomplete over the width of the test strip, or uneven in color, it should be interpreted as positive. A positive result indicates the presence of elevated MxA and/or CRP proteins.

Negative Result

If only a blue control line is visible in the **Result Window**, the test is deemed negative. A negative result indicates a lack of elevated MxA and CRP proteins.

Invalid Result

The absence of the blue control line indicates an invalid result.

If an invalid result occurs, the test must be discarded and the patient retested using a new FebrIDx test. Choose an alternative puncture site on a different finger when retesting the patient.

Notes:

- A blood fluid wave will migrate up the **Result Window** and gradually disappear as the test develops.
- Faint blood streaks may be visible along the sides of the **Result Window** and are acceptable for reading purposes.
- If the background of the **Result Window** has not cleared sufficiently for interpretation of results after 30 minutes, discard the test and retest the patient with a new FebrIDx test.
- FebrIDx test results are stable for up to one (1) hour. Do not interpret the test results after this period of time.

Control line

A blue control line must appear in the **Result Window** for the test to be valid. The absence of a blue control line indicates an invalid result.

External Controls

External controls may be used to demonstrate that the reagents and assay perform properly. FebrIDx external controls are available directly through Lumos Diagnostics. FebrIDx external controls consist of one (1) positive control and one (1) negative control.

For ordering external controls, refer to the Ordering and Contact Information section of this package insert. Refer to the FebrIDx external controls package insert for instructions on how to run the external controls. External controls will have an individual expiration date printed on each package. Discontinue use if external controls are expired.

LIMITATIONS

1. The FebrIDx test is best used within three (3) days of acute onset fever and seven (7) days of new onset respiratory symptoms.
2. Fresh capillary blood (fingerstick) must be used on the FebrIDx test.
3. Five (5) µl of blood must be applied to the test strip in order for the test to run properly. An erroneous result may occur if an insufficient blood sample is applied to the test.
4. The following conditions may lead to erroneous results:
 - Current immunosuppressive state or use of immunosuppressive drugs
 - Current use of oral anti-infective drugs
 - Current use of interferon therapy (e.g. for multiple sclerosis, HIV, hepatitis B/C)
 - Live viral immunization within the last 30 days
 - Major trauma, major surgical intervention, and severe burns within the preceding 30 days
 - Chronic fevers lasting more than 7 days
5. Many bacteria – including Streptococcus, Staphylococcus, Hemophilus, and enterobacter species – are known to colonize in the respiratory tract and their significance is unknown. Rhinovirus and coronavirus typically cause a non-invasive localized infection without an invasive systemic host response; while Herpes Simplex Virus, Epstein-Barr Virus, and Cytomegalovirus have a high frequency of DNA shedding. FebrIDx will not identify bacterial colonization, localized infections, or periodic viral shedding without an associated systemic host response.

EXPECTED VALUES

The prevalence of ARI varies during the year and from region to region, with outbreaks typically occurring during fall and winter. In developing countries, 30% of all patient consultations and 25% of all pediatric admissions are due to ARI.¹³ Additionally, fever has been shown to represent up to 30% of all pediatric office visits, making it the most common presenting sign of illness.^{14,15} ARIs are the leading cause of morbidity, accounting for: 20% of medical consultations, 30% of absenteeism, and 75% of all antibiotic prescriptions.¹⁶

PERFORMANCE CHARACTERISTICS

FebrIDx underwent two prospective, multicenter, blinded clinical trials with trained and untrained operators to determine the positive and negative agreement of the FebrIDx test at identifying a clinically significant immune response to diagnose and differentiate viral and/or bacterial community-acquired acute febrile upper respiratory infections as compared to the evaluation of expert clinical reviewers assessment of clinical standardized microbiologic and laboratory testing. The study was performed on patients older than 1 years that presented within 3 days of an acute onset fever and within 7 days of new onset respiratory symptoms consistent with a community-acquired acute febrile respiratory infection to primary care and urgent care outpatient offices and emergency departments.

Viral testing:

- Film Array polymerase chain reaction (PCR): Influenza A/B, Adenovirus, Respiratory Syncytial Virus 1-2, Parainfluenza virus 1-4, Metapneumovirus, Coronavirus, and Rhinovirus*
- Supplemental real-time PCR for Epstein-Barr Virus, Herpes Simplex Virus, and Cytomegalovirus**

Bacterial testing:

- Film Array PCR for atypical bacteria: *Chlamydomphila pneumoniae*, *Mycoplasma pneumoniae*, *Bordetella pertussis*, *Fusobacterium necrophorum*, *Neisseria gonorrhoeae*.***
- Oropharyngeal cultures (blood, chocolate, and MacConkey plates)

Laboratory testing:

- Procalcitonin (PCT) and white blood cell (WBC) count

Pharyngeal bacterial colonization was differentiated from true systemic bacterial infection if cell culture growth occurred in association with an elevated PCT level. Patients with microbiologically unconfirmed respiratory infection and a normal PCT level were deemed negative for clinically significant infection.

A reference-testing algorithm combined with an expert physician panel over-read served as the reference standard. Each patient underwent the following 6 diagnostic tests: (1) throat swab bacterial culture; (2) multiplex polymerase chain reaction (PCR) of a combined nasopharyngeal and oropharyngeal (NP/OP) sample using the Filmarray® Respiratory Panel (Biomérieux, Inc.; Marcy-L'Etoile, France); (3) real-time reverse transcriptase PCR of an NP/OP sample for Epstein-Barr virus (EBV); (4) EBV IgM serum antibody with the Immunosimplicity® IS-EBV-VCA IgM test kit (Diamedix Co; Miami Lakes, FL, USA); (5) serum PCT concentration measurement with the Brahms PCT Kryptor analyzer (Thermo Fisher; Waltham, MA, USA); (6) wbc count with band differential percentage, and (7) MxA protein ELISA and high sensitivity c-reactive protein enzyme immunoassay (Biocheck; Foster City, CA, USA). Reference testing was completed at a central laboratory and blinded to patients, treating clinicians, and study personnel who performed FebrIDx testing.

The reference testing algorithm classified patients as having a bacterial infection if any of the following 5 criteria were met: (1) throat culture positive for a bacteria that commonly causes pharyngitis (group A and C β-hemolytic Streptococci, N. gonorrhoeae, C. diphtheria, A. haemolyticum) plus PCT ≥0.1 ng/mL; (2) throat culture positive for any other bacteria plus PCT ≥0.15 ng/mL; (3) NP/OP sample PCR positive for atypical bacteria (M. pneumoniae, C. pneumoniae, B. pertussis) plus PCT ≥0.1 ng/mL; (4) PCT ≥0.25 ng/mL plus no identified pathogen; (5) PCT ≥0.15 ng/mL plus wbc ≥15,000 cells/mcl or presence of wbc bands plus no identified pathogen.

The reference testing algorithm classified patients as having a viral infection if any of the following 3 criteria were met: (1) NP/OP sample PCR positive for influenza A or B, Adenovirus, Respiratory Syncytial Virus (RSV), Human Metapneumovirus, Parainfluenza viruses 1-4,

QUALITY CONTROL

Procedural Controls

The FebrIDx test contains the following built-in procedural controls. For daily quality control, Lumos Diagnostics recommends documenting that these internal procedural controls are checked for the first sample tested each day.

Unused test

An unused FebrIDx test has faint orange lines in the **Result Window**, indicating the potential appearance of control and result lines.

Fluid wave

A blood fluid wave will migrate up the **Result Window** and gradually disappear as the test develops. Faint blood streaks may be visible along the sides of the **Result Window** and are acceptable for reading purposes.

Cytomegalovirus, and Herpes Simplex Virus; (2) NP/OP sample PCR positive for EBV plus serum IgM positive for EBV; (3) PCT between 0.15 ng/mL and 0.25 ng/mL plus wbc <15,000 cells/mcl plus no wbc bands plus no identified pathogen. Otherwise, the patients will be characterized as a viral infection. Patients who did not meet any of these criteria for bacterial or viral infection were classified as negative by the reference testing algorithm.

| Analyte Concentration | | Number Read as Positive | | Total tests | % Read as Positive | |
|-----------------------|----------|-------------------------|-----|-------------|--------------------|-----|
| MxA ng/mL | CRP mg/L | MxA | CRP | | MxA | CRP |
| 0 | 0 | 0 | 0 | 120 | 0 | 0 |
| 22 | 3 | 1 | 0 | 120 | 1 | 0 |
| 89 | 28 | 109 | 109 | 120 | 91 | 91 |
| 266 | 84 | 120 | 120 | 120 | 100 | 100 |
| 0 | 13 | 0 | 36 | 120 | 0 | 30 |
| 0 | 60 | 2 | 120 | 120 | 2 | 100 |
| 0 | 180 | 0 | 120 | 120 | 0 | 100 |
| 0 | 42 | 0 | 118 | 120 | 0 | 98 |

| Substance | Concentration |
|-------------------|-----------------|
| Ascorbic Acid | 175 mg/dL |
| Aspirin | 0.065 mg/mL |
| Conj. Bilirubin | 40 mg/dL |
| Dayquil® | 1:5000 Dilution |
| Free Bilirubin | 40 mg/dL |
| Glutathione | 300 µM |
| Hemoglobin | 20 g/dL |
| IL-6 | 1000 µg/L |
| IL-8 | 1000 µg/L |
| Interferon | 40 ng/mL |
| Isoniazid | 300 µM |
| L-Dopa | 300 µM |
| Oxaloacetate | 300 µM |
| Procalcitonin | 1000 µg/L |
| Rheumatoid Factor | 350 IU/mL |
| Robitussin® | 1:5000 Dilution |
| Triglycerides | 10 mg/mL |
| Tylenol® | 1:5000 Dilution |

PRECISION AND REPRODUCIBILITY STUDY

Samples were prepared in fresh EDTA whole blood with recombinant MxA and CRP proteins. Eight samples, consisting of a combination of C5 and C95 concentrations of the analytes namely, MxA and CRP were tested. A total of 960 determinations were performed by a combination of laboratory technicians, nurses, and clinicians at three different sites over five contiguous days during a two-week period. The study demonstrates overall reproducibility among three lots of material, among three separate sites, and among six separate users.

INTERFERING SUBSTANCES

The analytical specificity of the FebrIDx test was determined by evaluating a series of interfering substances mixed with samples with and without MxA and CRP antigens at their C95 (LoD) and C5 concentrations. Positive and negative interference of the following substances was examined. The list includes typical medications that febrile patients may be taking and the proteins and other substances normally found in blood.

The normal reference range for Rheumatoid Factor (RF) is less than 15 IU/mL. At a RF concentration of 350 IU/mL, positive interference in the presence of a C5 concentration was observed for MxA and CRP test lines. At a RF concentration of 150 IU/mL, positive interference in the presence of a C5 concentration was exhibited for CRP, but not for the MxA test line. No interference was observed for any test lines at a RF concentration of 116 IU/mL.

FEBRIDx PERFORMANCE DATA

Clinical Study N = 371 patients enrolled within 72 hours of a reported fever

- 165 Asymptomatic patients
- 205 Symptomatic URI patients
 - 38% confirmed infectious by gold standards
 - 12% Bacterial
 - 26% Viral
 - 62% Microbiologically Unconfirmed Respiratory Illness (MURI)
 - Two invalid tests occurred while testing asymptomatic patients
 - One URI patient had no cultures or blood work performed so this patient data was excluded
- **Bacterial infection**
 - Positive Agreement = 80% (20/25), 95% CI [58.7-92.4]
 - Negative Agreement = 93% (168/180), 95% CI [88.4-96.4]
 - Positive Predictive Value = 62% (20/32), 95% CI [43.7-78.3]
 - Negative Predictive Value = 97% (168/173), 95% CI [93.0-98.9]
- **Viral infection**
 - Positive Agreement = 86% (56/65), 95% CI [74.8-93.1]
 - Negative Agreement = 89%* (124/140), 95% CI [89.7-84.6]
 - Positive Predictive Value = 78% (56/72), 95% CI [66.2-86.4]
 - Negative Predictive Value = 93% (124/133), 95% CI [88.5-97.0]

*Bacterial colonization, localized infections, or periodic viral shedding without an associated systemic response was not detected. The presence of elevated PCT and/or WBC in association with culture growth was required to differentiate bacterial colonization from true systemic infection. Since rhinovirus and coronavirus are frequent colonizers of the respiratory tract and only cause a clinically significant infection in approximately 10% of patients, these two viruses were only considered true systemic infections if their presence was associated with elevated MxA; and therefore are not specifically called out in the intended use.

Rhinovirus was confirmed present by PCR in 52 subjects; however, only 8/52 patients demonstrated an elevation in MxA. Of those patients with confirmed rhinovirus and elevated MxA, FebrIDx correctly identified 5/8 patients with rhinovirus. Coronavirus was confirmed present by PCR in 6 subjects; however, only 1/6 patients demonstrated an elevation in MxA. FebrIDx correctly identified the 1 patient with confirmed coronavirus and elevated MxA.

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
REFERENCES

1. Nakabayashi M, Adachi Y, Itazawa T, et al. MxA-based recognition of viral illness in febrile children by a whole blood assay. *Pediatr Res* 2006;60(6):770-74.
2. Itazawa T, Adachi Y, Imamura H, et al. Increased lymphoid MxA expression in acute asthma exacerbation in children. *Allergy*. 2001;56(9):895-98.
3. Forster J, Schweizer M, Schumacher RF, et al. MxA protein in infants and children with respiratory tract infection. *Acta Paediatr* 1996;85(2):163-67.
4. Engelmann I, Dubos F, Lobert PE, et al. Diagnosis of viral infections using myxovirus resistance A (MxA). *Pediatrics* 2015;135:e985-93.
5. Halminen M, Ilonen J, Julkunen I, et al. Expression of MxA protein in blood lymphocytes discriminates between viral and bacterial infections in febrile children. *Pediatr Res* 1997;41(5):647-50.
6. Sambursky R, Shapiro N. Evaluation of a combined MxA and CRP point-of-care immunoassay to identify viral and/or bacterial immune response in patients with acute febrile respiratory infection. *Eur Clin Res J* 2015;10;2:28245.
7. Andreola B, Bressan S, Callegaro S, et al. Procalcitonin and c-reactive protein as diagnostic markers of severe bacterial infections in febrile infants and children in the emergency department. *Ped Infect Dis J* 2007;26(8):672-77.
8. Simon L, Gauvin F, Amre DK, et al. Serum procalcitonin and c-reactive protein levels as markers of bacterial infection: a systematic review and meta-analysis. *Clin Infect Dis*. 2004;39:206-17.
9. Falk G, Fahey T. C-reactive protein and community-acquired pneumonia in ambulatory care: Systematic review of diagnostic accuracy studies. *Fam Pract*. 2009;26(1):10-21.
10. Salonen EM, Vaheiri A. C-reactive protein in acute viral infections. *J Med Virol* 1981;8(3):161-7.
11. Smith RP, Lipworth BJ. C-reactive protein in simple community-acquired pneumonia. *Chest*. 1995;107(4):1028-31.
12. Smith RP, Lipworth BJ, Cree IA, et al. C-reactive protein. A clinical marker in community-acquired pneumonia. *Chest* 1995;108(5):1288-91.
13. Jain N, Lodha R and Kabra Sk. Upper respiratory tract infections. *Indian J Pediatr*. 2001;68:1135-1138.
14. Wright P, Thompson J, McKee K, et al. Patterns of illness in the highly febrile young child: epidemiologic, clinical and laboratory correlates. *Pediatrics* 1981;67:694-700.
15. Eskerud J, Laerum E, Fagerthun H, et al. Fever in general practice: I. Frequency and diagnoses. *Fam Pract* 1992;263-69.
16. Acute respiratory infections: the forgotten pandemic. Bulletin of the World Health Organization. 1998;76(1):101-107.

ORDERING AND CONTACT INFORMATION / INFORMACIÓN DE CONTACTO Y PARA PEDIDOS / BESTELINFORMATIE EN CONTACTGEGEVENS /
INFORMATIONS DE COMMANDE ET COORDONNÉES / BESTELL- UND KONTAKTINFORMATIONEN / INFORMAZIONI PER GLI ORDINI E CONTATTI

Manufacturer and United States Representative / Fabricante y representante en
Estados Unidos / Fabrikant en gemachtigde in de Verenigde Staten / Fabricant
et mandataire pour les États-Unis / Hersteller und Bevollmächtigter in den USA /
Produttore e rappresentante USA

European Representative / Representante europeo / Europees Gemachtigde /
Mandataire pour l'Europe / Bevollmächtigter in der EU / Rappresentante
per l'Europa

7227 Delainey Court, Sarasota, FL 34240 USA
t: +1.941.556.1850 f: +1.941.556.1851
RPSdetectors.com

Altenhofstr. 80, D-66386 St. Ingbert, Germany
t: +49.6894.581020 f: +49.6894.581021
Mt-procons.com

Lumos Diagnostics is a trade name of Lumos Diagnostics Holdings PTY LTD. Rapid Pathogen Screening, Inc. is a wholly owned subsidiary of Lumos Diagnostics Holdings PTY LTD. / Lumos Diagnostics es un nombre comercial de Lumos Diagnostics Holdings PTY LTD. Rapid Pathogen Screening, Inc. es una filial de titularidad exclusiva de Lumos Diagnostics Holdings PTY LTD. / Lumos Diagnostics is een handelsnaam van Lumos Diagnostics Holdings PTY LTD. Rapid Pathogen Screening, Inc. is een volle dochteronderneming van Lumos Diagnostics Holdings PTY LTD. / Lumos Diagnostics est une marque commerciale de Lumos Diagnostics Holdings PTY LTD. Rapid Pathogen Screening, Inc. est une filiale en propriété exclusive de Lumos Diagnostics Holdings PTY LTD. / Lumos Diagnostics ist eine Handelsbezeichnung von Lumos Diagnostics Holdings PTY LTD. Rapid Pathogen Screening, Inc. ist ein Tochterunternehmen von Lumos Diagnostics Holdings PTY LTD. / Lumos Diagnostics è un nome commerciale di Lumos Diagnostics Holdings PTY LTD. Rapid Pathogen Screening, Inc. è una consociata interamente controllata da Lumos Diagnostics Holdings PTY LTD.

©2020 Lumos Diagnostics. All rights reserved / Reservados todos los derechos. / Alle rechten voorbehouden / Tous droits réservés. / Alle Rechte vorbehalten / Tutti i diritti riservati