## INSTRUCTIONS FOR USE

# Actim<sup>®</sup> PROM

EN	Instructions for use	3
CS	Návod k použití	6
DA	Brugervejledning	9
DE	Gebrauchsanweisung	12
ES	Instrucciones de uso	15
FI	Käyttöohje	18
FR	Notice d'utilisation	21
IT	Istruzioni d'uso	24
NL	Gebruiksaanwijzing	27
NO	Bruksanvisning	30
PL	Instrukcja użytkowania	33
PT	Instruções de utilização	36
RO	Instrucțiuni de utilizare	39
RU	Инструкция по применению	42
SV	Bruksanvisning	45







# Structure of dipstick

[1] Dip area

[2] Result area

## Test procedure and results

# Instructions for use

Numbers [1] - [10] refer to illustrations on inner cover.

## Structure of dipstick

[1] Dip area [2] Result area [3] Test line [4] Control line

## Intended use

Actim<sup>®</sup> PROM test is a visually interpreted, qualitative immunochromatographic dipstick test for detection of amniotic fluid in vaginal secretions during pregnancy. Actim PROM detects IGFBP-1, which is a major protein in amniotic fluid and a marker of the amniotic fluid in a vaginal sample. The test is intended for professional use to help diagnose the rupture of fetal membranes (ROM) in pregnant women.

## Kit components

Actim PROM kit 30831ETAC contains 10 and the kit 30832ETAC contains 20 test packs with instructions for use. The kit 30830ETAC contains one test pack and instructions for use.

The components of each Actim PROM test pack (30821ETAC) are:

- One sterile polyester swab for specimen collection.
- One tube of Specimen Extraction Solution (0.5 ml). This phosphate-buffered solution contains bovine serum albumin (BSA), protease inhibitors and preservatives.
- One dipstick in a sealed aluminum foil pouch with desiccant.

## Storage

Store the test kit at +2...+25 °C. Stored unopened, each component can be used until the expiry date marked on the component. The kit can also be stored for 2 months at +2...+30 °C. Use the dipsticks shortly after their removal from the aluminum foil pouch.

## Specimen collection

The specimen is vaginal secretion that is extracted into the Specimen Extraction Solution provided. A sample is obtained using a sterile polyester swab (provided in the kit). The sample should be collected prior to performing digital examination and/or transvaginal ultrasound. Take care not to touch anything with the swab before taking the sample. Separate the labia and carefully insert

the tip of the swab into the vagina toward the posterior fornix until resistance is met. Alternatively the sample can be taken from the posterior fornix during a sterile speculum examination. The swab should be left in the vagina for 10–15 seconds to allow it to absorb the vaginal secretion [5].

Open the Specimen Extraction Solution tube and put it in a vertical position. The specimen is extracted immediately from the swab by swirling the swab vigorously in the extraction solution for 10–15 seconds [6]. Press the swab against the wall of the Specimen Extraction Solution tube to remove any remaining liquid from the swab. Discard the swab after extraction. **NOTE!** Do not leave the swab to the tube.

Specimens should be tested as soon as possible after extraction but in any case no more than 4 hours after specimen collection and extraction. If a specimen cannot be tested within this time it should be frozen. After thawing, the specimens should be mixed and tested as described below.

## Test procedure and interpretation of the results

1. If stored refrigerated, allow the aluminum foil pouch and the Extraction Solution tube to reach room temperature. Open the foil pouch containing the dipstick by tearing. Do not touch the yellow dip area at the lower part of the dipstick. Identifying marks may be written on the upper turquoise part of the dipstick. The dipstick must be used shortly after its removal from the foil pouch.
2. Place the yellow dip area into the extracted sample [7] and hold it there **until you see the liquid front enter the result area [8]**. Remove the dipstick from the solution and place it in a horizontal position.
3. The result can be interpreted as positive as soon as two blue lines become visible in the result area. Negative result should be read at **5 minutes [9]**. **Do not pay attention to any lines appearing later than 5 minutes.**
4. If **two blue lines, the test line and the control line**, appear, the test result is **positive**. If **one blue line, the control line**, appears, the test result is **negative**.

If **the control line does not appear**, the test is **invalid [10]**.

## Limitations of the test

- The test is intended for *in vitro* diagnostic use only.
- If rupture of fetal membranes has occurred **but the leakage of amniotic fluid has ceased** more than 12 hours before the specimen is taken, IGFBP-1 may have been degraded by proteases in the vagina and the test may give a negative result.
- A negative test result is an indication on the present condition and cannot be used to predict the forthcoming.

## Notes

- No quantitative interpretation should be made based on the test results.
- The test requires about 150 µl of extracted sample to ensure proper performance of the test.
- Care must be taken when placing the dipstick in the sample tube. The upper part of the dipstick must stay dry.
- Do not use a dipstick that has become wet before use, because moisture damages the dipstick.
- Do not use a dipstick if you notice a blue coloring in the result area before testing.
- Do not use the dipstick if its aluminum foil pouch or the seals of the pouch are not intact.
- Use only the swab provided with the kit.
- Improper sampling may lead to false negative result.
- When dipping, be careful to hold the dipstick in position (with the dip area in the sample extract) until the sample liquid front reaches the result area.
- If the control line does not appear, the test is invalid, and should be repeated using another dipstick.
- If the test result cannot be interpreted clearly it is recommended that the test be repeated.
- At five minutes the appearance of any faint-to-dark blue test line along with a control line indicates a positive result. However, do not pay attention to any lines appearing after 5 minutes.

- If only the control line is visible, the result should be interpreted as negative only after 5 minutes have elapsed.
- A positive Actim PROM test result, although detecting the presence of amniotic fluid in the sample, does not locate the site of the rupture.
- As with all diagnostic tests, results must be interpreted in the light of other clinical findings.
- All biological specimens and materials must be treated as potentially hazardous, and disposed of in accordance with local authority guidelines.

## Principle of the test

The concentration of IGFBP-1 (insulin-like growth factor binding protein-1) in amniotic fluid is 100 to 1000 times higher than in maternal serum. IGFBP-1 is not usually present in the vagina, but after rupture of fetal membranes, amniotic fluid with a high concentration of IGFBP-1 mixes with vaginal secretions. In Actim PROM test, a specimen of vaginal secretion is taken with a sterile polyester swab and the specimen is extracted into Specimen Extraction Solution. The presence of IGFBP-1 in the solution is detected using a dipstick.

The test is based on immunochromatography. It involves two monoclonal antibodies to human IGFBP-1. One is bound to blue latex particles (the detecting label). The other is immobilized on a carrier membrane to catch the complex of antigen and latex-labeled antibody and indicate a positive result. When the dip area of the dipstick is placed in an extracted sample, the dipstick absorbs liquid, which starts to flow up the dipstick. If the sample contains IGFBP-1 it binds to the antibody labeled with latex particles. The particles are carried by the liquid flow and, if IGFBP-1 is bound to them, they bind to the catching antibody. A blue line (test line) will appear in the result area if the concentration of IGFBP-1 in the sample exceeds the detection limit of the test. A second blue line, the control line, confirms correct performance of the test.

### Performance of the test

The lowest detectable amount of IGFBP-1 in the extracted sample is 25 µg/l (calibrated against Behring Institut IGFBP-1 [PP12] preparation Lot 307/323).

### Analytical Sensitivity

The analytical sensitivity (detection limit) of the Actim PROM test was identified by evaluating different concentrations of IGFBP-1 in extracted samples on three different lots of the Actim PROM test. Two different operators each interpreted ten devices run at each concentration under various lighting conditions for a total of 60 determinations per level. The Actim PROM test limit of detection is approximately 25 µg/l in extracted sample. The measuring range of the Actim PROM is approximately 25-500000 µg/l in extracted sample.

### Analytical Specificity

Analytical specificity (cross-reactivity) was tested with human IGFBP proteins at concentrations ranging from 10-5000 µg/l of each protein in extracted sample using three different lots of the Actim PROM test. No cross-reactivity was seen using human IGFBP-2, -3, -4, -5 and -6 proteins. The Actim PROM test is specific to human IGFBP-1.

### Repeatability and reproducibility

A panel of specimens consisting of samples of different IGFBP-1 concentration levels was evaluated for intra- (repeatability) and inter- (reproducibility) assay precision. For repeatability the samples were tested with 10 replicates during the same day using three different lots of the Actim PROM test. Repeatable results were obtained. For reproducibility a study of the Actim PROM test was conducted at three separate sites. Test operators (n=9) tested each level on five different days using one lot of the Actim PROM test. A total of 360 tests were performed (120 per site) with a total of 45 tests per sample type. The overall reproducibility of the Actim PROM test is 97 % [350/360] with no significant differences within runs (replicates tested by one operator), between runs (five different days), between sites (three sites) or between operators (nine operators).

### Diagnostic performance

Actim PROM test has been evaluated in several clinical studies. Selection of typical studies is shown in FIG 1 on the inner back cover.

### Interference testing

The following substances, conditions and micro-organisms were tested with Actim PROM test and were found not to affect Actim PROM test performance.

TESTED SUBSTANCE/CONDITION/ MICRO-ORGANISM	INTERFERENCE (NO/YES)
Whole blood	No
Semen	No
Urine	No
pH 3.5-8.5	No
Ultrasound Transmission Gel	No
Pevaryl (active ingredient: econazol.nitras)	No
Gyno-Trosyd (tioconazol)	No
Flagyl (metronidazole)	No
Canesten (clotrimazol)	No
Personal Lubricant	No
Baby Oil	No
Baby Powder	No
Feminine Deodorant Suppositories	No
Vaginal Gel	No
Shower and bath products	No
Candida albicans	No
Gardnerella vaginalis	No
Neisseria gonorrhoea	No
Chlamydia trachomatis	No
HSV-1	No
HSV-2	No

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**Erdemoglu and Mungan T.**

Significance of detecting insulinlike growth factor binding protein-1 in cervicovaginal secretions: comparison with nitrazine test and amniotic fluid volume assessment. Acta Obstet Gynecol Scand [2004] 83:622-626.

**Jain K and Morris PG.**

A clinical study to evaluate the usefulness of the MAST test in diagnosing pre-labour rupture of membranes. J Obstet Gynaecol [1998] 18:33-36.

**Kubota T and Takeuchi H.**

Evaluation of insulin-like growth factor binding protein-1 as a diagnostic tool for rupture of the membranes. J Obstet Gynecol Res [1998] 24:411-417.

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**Rutanen E-M et al.**

Evaluation of a rapid strip test for insulinlike growth factor binding protein-1 in the diagnosis of ruptured fetal membranes. Clin Chim Acta [1996] 253:91-101.

**Rutanen E-M, Pekonen F, Kärkkäinen T.**

Measurement of insulin-like growth factor binding protein-1 in cervical/vaginal secretions: comparison with the RO Mcheck Membrane Immunoassay in the diagnosis of ruptured fetal membranes. Clin Chim Acta [1993] 214:73-81.

# Actim<sup>®</sup> PROM

Actim Oy  
Klovinpellontie 3, FI-02180 Espoo, Finland  
Tel. +358 9 547 680  
[www.actimtest.com](http://www.actimtest.com)  
[actim@actimtest.com](mailto:actim@actimtest.com)

© Actim Oy. All rights reserved.  
Actim is a trademark of Actim Oy.

AOACE30831-2 2020-07-07