

# INTRODUCING DRIVECHECK


*Why take a chance?*

SEMI-QUANTITATIVE DISPOSABLE BREATHALZYER

DriveCheck is a single-use, semi- quantitative breath alcohol test designed to rapidly detect the presence of alcohol in exhaled breath, with relative BrAC levels of 0.00%, 0.02%, 0.05% and 0.08%.

## Benefits:

- Rapid results in just minutes.
- Detection at 4 levels (semi-quantitative)
- Easy to interpret
- Low cost point
- Discrete and ideal for any setting



DRIVECHECK™

  
SAFETYCHECK™  
Have fun, keep safe.

  
SPECTRUM  
MEDICAL DIAGNOSTICS