

ALLTEST Rheumatoid Factor Rapid Test

For Reliable Diagnosis of Rheumatoid Arthritis

- ◆ **Rheumatoid arthritis (RA)** is a long-term autoimmune disorder that primarily affects joints. Onset is most frequent during middle age and women are affected 2.5 times as frequently as men.
- ◆ RA typically results in warm, swollen, and painful joints. It is diagnosed by reviewing symptoms, conducting a physical examination, doing X-rays and lab tests.
- ◆ **Rheumatoid factor (RF)** is an autoantibody most relevant in RA. It is part of the usual disease criteria of RA. High levels of rheumatoid factor (generally above 20 IU/ml) occur in RA (present in 80%).
- ◆ ALLTEST™ Rheumatoid Factor (RF) Rapid Test Cassette utilizes RF antibodies to specifically detect the presence of RF in human serum/plasma without the use of an instrument.



ALLTEST Rheumatoid Factor Rapid Test Offers :

✓ High Accuracy

Sensitivity: 98.0%

Specificity: 97.7%


Accuracy: 97.8%

✓ Fast Results in 10 Minutes

✓ Simple Procedure



✓ Convenient Procedure



✓ Ordering Information

Cat. No.	Product Description	Specimen	Format	Pack	Cut-off	CE Status
ORF-302	Rheumatoid Factor Rapid Test	Serum/Plasma	Cassette	20 T	25IU/mL	CE

Hangzhou AllTest Biotech Co.,Ltd.