

PSA Rapid Testing



Patient-Centred Care
Digital Innovation

SpecCheck Rapid PSA Test (with Results+ Patient Reporting)

The PSA Rapid Test Cassette (Whole Blood/Serum/Plasma) is a professional in vitro diagnostic test for the qualitative detection of Prostate Specific Antigen (PSA).

Our new Results+ platform has been integrated with this test to quickly and efficiently provide test results to your consumers on the spot!



Identify whole blood specimens with 99.1% Accuracy.



Obtain results in just 10 minutes.



Determine PSA levels from >3ng/mL to <10ng/mL.

Patient-Centric Services

Improve patient care by offering value added screening services in-store. Pharmacists can help address current backlogs within hospitals & clinics and minimize the strain on the Canadian healthcare system

Digital Innovation with Results+

As healthcare delivery evolves and digital innovation persists, we want to help bridge the gap between medical device and software. Capture test images and print and/or email results with one click of a button!

Accuracy & Reliability

Excellent correlation to laboratory results. Quick, accurate & reliable results for enhanced healthcare.