



## **Time-resolved Fluorescence Immunoassay Analyzer (OG-G200)**

Test Markers for diabetes, inflammation, hormones, cardiac markers and more.



**Accurate results  
within eight  
seconds.**

## **The Benefits of Fluorescence Immunoassay Testing**

### **Accessibility & Convenience**

Easily accessible FIA testing in pharmacies, especially beneficial in rural and underserved areas.



### **Enhanced Healthcare Efficiency**

Pharmacies offering FIA testing reduce costs and hospital workload, boosting care efficiency during critical times.


## Ultra-Sensitive Detection

Capable of identifying minute biomarkers for brain injuries as low as 1.3 ng/mL—akin to finding a pinch of salt in a swimming pool.



## High Precision Cardiac Diagnostics

Fluorescent Immunoassay technology demonstrates less than 2.5% variation in results, ensuring consistent outcomes for accurate heart damage assessments.




## Secure Patient Wellness:

Bypass delays, embrace rapid health solutions.

## Swift Diagnosis, Immediate Relief:

Skip the lab wait, start treatment today.



## Enhance Healthcare Efficiency:

Rapid testing in emergency departments has shown a **65.6% reduction** in patient isolation time, significantly boosting treatment capacity and reducing resource use.

## Advance Healthcare Outcomes:

Enhance Precision and efficiency with Fluorescence Immunoassay Screening!

[Learn More](#)



### Specifications


Power supply	12V, 5A
Dimensions	232mm * 269mm * 147 mm
Display	2050 g
Weight	7 inch LCD Capacitive Screen, Res: 1024 * 600
Printer	Built-in thermal printer
Network	Supports Wi-Fi connection
Detection speed	8 seconds
Storage Capacity	10,000 sets of data
Language Interface	Chinese, English, other languages can be expanded


*Suitable for Clinics,  
Laboratories, Hospitals  
and Pharmacies.*



Ailment	Popular Test Cartridges Available
Cardiovascular Diseases	cTnl, CK-MB, Myoglobin Tests
Diabetes	HbA1c Test
Respiratory Diseases	Influenza A&B Test, RSV Test
Infectious Diseases	SARS-COV-2 & Influenza A/B Combo Test
Cancer	AFP, CEA Tests
Chronic Kidney Disease	Microalbumin, $\beta$ 2-Microglobulin Tests
Hypertension	Angiotensin I & II Tests
Thyroid Disorders	TSH, Free T3, Free T4 Tests
Arthritis	CRP, SAA Tests
Osteoporosis	$\beta$ -CTX, PINP, Osteocalcin Tests

For more information on our products and how they can help manage health effectively, please contact us.

 [info@spectrummdx.com](mailto:info@spectrummdx.com)

 (866) 287-2425

 <https://spectrummdx.com>